



Identifying Select Agents

Shannon Sharp
Bioterrorism Training Coordinator
MDHHS Bureau of Laboratories



Objectives

1. Define the responsibilities of a sentinel laboratory.
2. Describe the types of samples to collect to rule out select agents.
3. Determine the tests to perform to rule out suspected select agents.
4. List the steps to take once a select agent cannot be ruled out.




Laboratory Response Network for Biological Threats (LRN-B)

- ▶ The LRN-B was founded in 1999 by the Centers for Disease Control and Prevention (CDC), Federal Bureau of Investigation (FBI) and the Association of Public Health Laboratories (APHL) to coordinate laboratory response to biological, chemical, radiological threats and other high priority public health emergencies, including emerging infectious diseases



LRN-B cont.

- ▶ **National Laboratories**, including the CDC, U.S. Army Medical Research Institute of Infectious Diseases (USAMRIID), and the Naval Medical Research Center (NMRC), are responsible for specialized strain characterization, bioforensics, select agent activity and handling of highly infectious biological agents.
- ▶ **Reference Laboratories**, are responsible for investigation and confirmatory testing.
- ▶ **Sentinel Laboratories**, comprised of hospital-based and commercial laboratories, are responsible for the early detection and the rule-out or referral of potential biothreat agents.





Sentinel Laboratory Definition

The laboratory is certified to perform high complexity testing under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) by the Centers for Medicare & Medicaid Services (CMS) for the applicable Microbiology specialty or the laboratory is a Department of Defense (DoD) Laboratory certified under the DoD Clinical Laboratory Improvement Program or the laboratory is a veterinary medical diagnostic laboratory that is fully accredited by the American Association of Veterinary Laboratory Diagnosticians (AAVLD).

Laboratory in-house testing includes Gram stains and at least one of the following: lower respiratory tract, wound or blood cultures.*

* American Society for Microbiology (asm.org)



Responsibilities of Sentinel Lab

- Familiar with reportable disease guidelines in its jurisdiction; has policies and procedures to refer specimens or isolates suspected to contain biothreat agents to the local/state public health laboratory.
- Ensures personnel meet applicable federal regulations for packaging and shipping of infectious substances.
- Has policies and procedures for referral of suspect biothreat agent specimens and/or isolates reflecting the American Society for Microbiology (ASM) Sentinel Level Clinical Laboratory Protocols for Suspected Biological Threat Agents and Emerging Infectious Diseases.



Responsibilities of Sentinel Lab cont.

- Maintains capability to perform testing outlined in the ASM Sentinel Level Clinical Laboratory Protocols for Suspected Biological Threat Agents and Emerging Infectious Diseases and demonstrates annual competency by participation in proficiency testing or exercises.
- Has a Class II or higher Certified Biological Safety Cabinet (BSC).
- Complies with Biological Safety Level II (BSL-2) practices as outlined in the current edition of the Biosafety in Microbiological and Biomedical Laboratories guidelines (BMBL).

Responsibilities of Sentinel Lab cont.

- Complies with applicable Occupational Safety and Health Administration (OSHA) regulations for a respiratory protection program.
- Complies with the rules and regulations of the Federal Select Agent Program.





American Society for Microbiology (ASM)

<https://www.asm.org/>

Select Agent Guidelines Available on the ASM Website

Anthrax (<i>Bacillus anthracis</i>)	Novel Influenza Viruses
<i>Brucella</i>	Plague (<i>Yersinia pestis</i>)
Botulinum Toxin	Smallpox
<i>Burkholderia</i>	Staphylococcal Enterotoxin B
<i>Coxiella burnetii</i>	Tularemia (<i>Francisella tularensis</i>)

Additional guidelines: BT Readiness Plan & Packaging and Shipping



Select Agents: Information & Identification


Bacillus anthracis “Anthrax”

- Zoonotic organism (primarily herbivores)
- Soil reservoir
 - Also found in water and plants
- No person to person transmission
- Safety Consideration:
 - As soon as *B. anthracis* is suspected in the laboratory, perform ALL further work within a Class II Biological Safety Cabinet (BSC) using BSL-3 practices.



Bacillus anthracis

- Infection Types–Based on Entry Into Body:
 - Inhalation (most deadly form of anthrax)
 - Inhale spores which can occur following intentional aerosol release
 - I.e. 2001 anthrax biocrime
 - If untreated, 85 to 90% patients will not survive
 - Cutaneous
 - Most common form of naturally occurring cases
 - Spores enter skin through cuts and abrasions
 - Gastrointestinal
 - Consume contaminated undercooked meat
 - Injection
 - Identified in heroin–injecting drug users in northern Europe.






Sample Collection for *Bacillus anthracis*

A. Collection and Transport of Clinical Specimens for Laboratory Rule-Out Testing*		
Cutaneous	Vesicular (early) stage	Unroof vesicle and aspirate fluid or collect with two sterile swabs (dacron)
	Eschar (late) stage	Insert swab (dacron) beneath the edge of the eschar, rotate swab or obtain an aspirate. Transport specimens at room temperature.
Gastrointestinal		Stool (> 5 grams...pecan size), collect and transport in a leak proof sealed container
		Collect blood (late stage of infection) directly into an appropriate blood culture bottle (aerobic and anaerobic)
		Transport specimens and bottles at room temperature
Inhalational		Sputum
		Blood: collect directly into an appropriate blood culture bottle (aerobic and anaerobic)
		Cerebral Spinal Fluid only if signs of meningitis occur
		Transport specimens and bottles at room temperature
Postmortem Tissue		Tissue pieces should be collected and kept moist
		Transport in sterile container at room temperature within 1 hour of collection

Bacillus anthracis

- Gram Stain
 - Large Gram positive rods
 - Spores may be found in cultures grown in 5% CO₂ but are not usually seen in clinical samples

Gram stain of *Bacillus anthracis*



Gram positive, endospore-forming
Bacillus anthracis



Bacillus anthracis

▶ Colony Morphology

- Grows well on Blood (BAP) & Chocolate (CHOC) agar
 - Non-hemolytic on BAP
- Does not grow on MacConkey (MAC)

24 HOURS ON BAP @ 35° C



Colonies on BAP showing
"Medusa Head" morphology



Bacillus anthracis

➤ Additional Tests

- Catalase: Positive
- Motility: Non-motile

➤ Common Misidentifications

- May not be identified in common automated ID systems, including MALDI TOF, and possible misidentifications may include: *Bacillus megaterium* and other *Bacillus* species.


Note: *Bacillus cereus* Group includes *B. anthracis*, but automated ID systems may not alert microbiologist beyond this group identification.

SAFETY: As soon as *B. anthracis* is suspected, perform ALL further work in a Class II BSC using BSL-3 practices.

Major characteristics of *Bacillus anthracis*:

Gram stain morphology: Large, Gram positive rods. Spores may be found in cultures grown in 5% CO₂ but not usually in clinical samples

Colony morphology: Ground glass appearance, non-pigmented, gamma hemolytic (no hemolysis) on BAP,
No growth on MAC (or EMB)



Bacillus cereus biovar *anthracis*

- Centers for Disease Control and Prevention (CDC) added organism to the Tier 1 select agent list on Oct. 14, 2016
- First described in gorillas and chimps in Cameroon & Côte d'Ivoire
 - Has been isolated from elephants & goats in Africa
 - No human infections to date
- Genetically similar to *B. anthracis* and produces all of the primary *B. anthracis* virulence factors

Bacillus cereus biovar *anthracis*

Characteristic	<i>B. anthracis</i>	<i>B. cereus</i>	<i>B. cereus</i> biovar <i>anthracis</i> CI ¹	<i>B. cereus</i> biovar <i>anthracis</i> CA ²
Hemolysis ³	-	+	-	-
Motility ⁴	-	+	+/-	+/-
Gamma phage susceptibility ⁵	+	-	-	-
Penicillin G ⁶	S	R	S	R
Capsule	+	Absent in vitro	+	+

▶1:CI = Côte d'Ivoire strains,
from chimpanzees

▶2:CA = Cameroon strains
from gorillas/chimpanzees

▶3:Hemolysis:
◦+=beta hemolytic on SBA; -=non-hemolytic


▶4: Motility:+ = motile; - = non-motile

▶*B. cereus* biovar *anthracis* strains are usually motile with exception of those recovered from goat strains from Democratic Republic of the Congo which were non-motile (3)

▶5: Gamma phage susceptibility:
▶ += susceptible; - = resistant

▶6. S=susceptible; R=resistant

24 h growth on BAP, 5% CO₂



Cameroon (CA) strain



Côte d'Ivoire (CI) strain

Bacillus cereus biovar *anthracis*

- This organism is very rare and travel history plays a key role if the physician suspects this organism.
- Recommendations:
 - Follow local public health guidelines to assess whether the public health lab or clinical lab should contact the patient's physician to determine likely clinical significance
 - Does the patient have an anthrax-like clinical syndrome
- If *Bacillus cereus* biovar *anthracis* is suspected, contact your LRN laboratory for further guidance.

Brucella spp.

- Found worldwide
 - Especially found in countries that lack effective public health & domestic animal health programs
 - I.e. Mediterranean Basin, Mexico, Eastern Europe
- Modes of transmission:
 - Ingestion—Most common method
 - Eating unpasteurized dairy products
 - Inhalation
 - Laboratorians—**Most commonly reported laboratory-associated bacterial infection**
 - Cutaneous & Mucous Membrane
 - Contact with infected animals



Brucella spp.

Red arrows indicate select agent species

Host Animals for <i>Brucella</i> Species	
Species	Main Animal Host(s)
<i>B. abortus</i>	Cattle http://www.cdc.gov/brucellosis/veterinarians/cattle.html
<i>B. canis</i>	Dogs http://www.cdc.gov/brucellosis/veterinarians/dogs.html
<i>B. ceti</i>	Cetaceans (dolphins, porpoises, whales)
<i>B. melitensis</i>	Goats, sheep, camels
<i>B. microti</i>	Common vole
<i>B. neotomae</i>	Wood rats
<i>B. ovis</i>	Sheep, goats
<i>B. pinnipediae</i>	Pinnipeds (seals, sea lions, walruses)
<i>B. suis</i>	Pigs

Brucella spp.

- Safety considerations
 - As soon as *Brucella* is suspected in the laboratory, perform ALL further work within containment such as a Class II Biological Safety Cabinet (BSC) and follow BSL-3 practices.



Brucella spp.

- A number of factors contribute to the risk of an accidental *Brucella* laboratory exposure, including:
 - Lack of experience working with the organism
 - Unknown or unidentified samples that arrive for analysis
 - Work performed on a *Brucella* isolate on an open bench, not under containment such as within a Class II BSC, or in BSL-3 conditions.

Brucella spp.



- **Minimize Exposure Risks:**
- Manipulate isolates of small gram-negative or gram-variable coccobacilli /rods within a Biological Safety Cabinet (BSC).
- Use primary barriers such as sealed safety centrifuge cups, PPE, and Class II or higher BSC for procedures with a high likelihood of producing aerosols.
- Restrict access to the laboratory when work is being performed
- Minimize the creation of splashes or aerosols.
- Prohibit sniffing of opened culture plates

Brucella spp.

Collection and Transport of Clinical Specimens for Laboratory Rule-Out Testing	
Bone marrow or whole blood	Considered the best specimen for culture but can be limited in chronic cases. Collect directly into an appropriate blood culture bottle. Aseptically inoculate liquid blood culture bottles with maximum amount of blood or body fluid per manufacturers' instructions. Incubate at 35°C. Transport bottles at room temperature as soon as possible to obtain the diagnosis. Alternatively, follow the manufacturer's instructions for the lysis-centrifugation method and inoculate pellet to BAP, CHOC and MAC. Incubate plates at 35°C in a humidified incubator with 5 to 10% CO ₂ .
Joint or abdominal fluid	Collect directly into an appropriate blood culture bottle. Transport bottles at room temperature as soon as possible to obtain the diagnosis.
Spleen, liver abscesses	Tissue pieces (at least the size of a pea) should be collected and kept moist. Transport in sterile container at room temperature within 1 hour of collection. May add 1-2 drop of sterile saline to keep moist.
Serum	Collect at least 1 mL without anticoagulant for serologic diagnosis. ¹ Store at 4°C until testing is performed. Acute specimen is collected as soon as possible after onset of disease. Convalescent-phase should be collected >14 (14 – 21) days after the acute specimen.
Tissues	Inoculate BAP, CHOC and MAC and incubate at 35°C in a humidified incubator with 5 to 10% CO ₂ for 7 days (MAC need only be incubated for 3 days at 35°C in ambient air or 5 to 10% CO ₂ .) Humidity may be maintained by placing a pan of water in the bottom of the incubator or by wrapping the plates with gas permeable tape.

Brucella spp.

- ▶ Gram Stain
 - Small gram negative coccobacilli or tiny GNR
 - Sometimes will take up the crystal violet and appear as gram positive




Brucella spp.

➤ Colony Morphology

- Very slow growing (may look like haze at 24 h)
- Can take 2–4 days to grow in blood cultures
- Grows on BAP or CHOC; No growth on MAC

24 Hours on BAP



48 Hours on BAP




Brucella spp.

➤ Biochemical Tests


- Catalase positive
- Oxidase positive
- Urea positive
- Satellite test negative

➤ Sometimes misidentification include:

- *Haemophilus spp.*
- *Oligella ureolytica*
- *Psychrobacter phenylpyruvicus*
- *Psychrobacter immobilis*
- *Bordetella bronchiseptica*



Brucella growing on the entire plate.
Negative for X and V factors.



Haemophilus demonstrating
growth only around the
Staphylococcus aureus streak.

Major characteristics of *Brucella* species:

Gram stain morphology: Small (0.4 x 0.8µm), Gram negative coccobacillus

THINK BRUCELLA

Colony morphology: Subculture positive aerobic blood culture to BAP, CHOC. Incubate in 5-10% CO₂ at 35°C, Spot BAP with *S. aureus* ATCC 25923 for satellite test. Note poorly growing colonies after 24 hour incubation on BAP and CHOC.

Incubate plates for at least two additional days if no growth in 24 hours.

Organism does NOT grow on MAC/EMB. (Pinpoint colonies infrequently after extended Incubation; 7 days)

Is the organism only growing on BAP without the need to satellite around the *S. aureus* at 24-48 hours?

No

Consider *Haemophilus*

Yes

Oxidase positive and catalase positive?

No

Think *Francisella* (see procedure)

Yes

Urea positive?

No

Reincubate and see written procedure

Yes

Brucella not ruled out.

Contact LRN Reference Level Laboratory to refer the isolate.

Report: Possible *Brucella* species submitted to LRN Reference Laboratory.

Antimicrobial therapy: Rifampin or Streptomycin plus Doxycycline

Burkholderia mallei (Glanders)

- Equines are primary reservoir but can effect other animals
- Used during WWI as biowarfare agent
- 4 types of infection:
 - Localized
 - Pulmonary
 - Bloodstream
 - Chronic




Burkholderia mallei

Collection and Transport of Clinical Specimens for Laboratory Rule-Out Testing	
Bone marrow or whole blood	<ul style="list-style-type: none"> • Considered the best specimen for culture • Collect directly into an appropriate blood culture bottle • Transport bottles at room temperature as soon as possible to obtain the diagnosis
Sputum or bronchoscopically obtained specimens	<ul style="list-style-type: none"> • Collect expectorated specimen into sterile transport cup or collect during bronchoscopy procedure. • Transport at room temperature up to 2h • If it is known that material will be transported from 2-24 h after collection, then store and transport at 2-8°C.
Tissue specimens (biopsies, abscess aspirates) and wound swabs	<ul style="list-style-type: none"> • Tissue pieces (at least the size of a pea) should be collected and kept moist • Transport in sterile container at room temperature within 1 hour of collection • Alternatively a swab from a tissue sample can be submitted in hospital transport tube with medium to stabilize specimen (e.g. Amies charcoal).
Urine	<ul style="list-style-type: none"> • Collect at least 1 mL into leak-proof container • Transport at room temperature up to 2 h • Refrigerate 2 up to 24 h until culture inoculation

Burkholderia mallei

➤ Gram Stain

- Small straight or slightly curved Gram negative coccobacilli with rounded ends
- Cells arranged in pairs, parallel bundles, or the Chinese letter form



Burkholderia mallei

➤ Colony Morphology

- On BAP, pinpoint to small grey colonies at 24h that may become smooth, grey, and translucent at 48h with no distinctive odor
- Non-hemolytic on BAP
- No growth or pinpoint colorless colonies (which may uptake light pink dye) on MAC after 48h

24 HOURS ON BAP




48 HOURS ON BAP



Burkholderia mallei

➤ Biochemical Tests

- Catalase positive
- Oxidase variable; most are negative
- Spot indole negative
- Non-motile (*Recommend tube test not wet mount due to potential aerosol production*)
- No growth at 42 °C
- Polymyxin B and colistin no zone, penicillin resistant, amoxicillin-clavulanate susceptible



Burkholderia mallei

Possible Misidentifications for *Burkholderia mallei* include:


Organism	Differential Test
<i>Burkholderia cepacia</i>	Resistant to amoxicillin-clavulanic acid, lactose fermenter (LF) on MacConkey and EMB, motile, arginine negative
<i>Chromobacterium violaceum</i>	Hemolysis, violet pigment on BAP, motile
<i>Pseudomonas stutzeri</i>	Growth on MacConkey, arginine negative
<i>S. maltophilia</i>	Growth on MacConkey, arginine negative
<i>Bacillus</i> spp. may appear Gram negative	Sensitive to penicillin
<i>Pandoraea</i> spp.	Growth on MacConkey
<i>Ralstonia</i> spp.	Growth on MacConkey

Major characteristics of *Burkholderia mallei*:

Gram stain morphology: Gram negative coccobacilli or small rods

Colony morphology: Poor growth at 24 hr; better growth of gray, translucent colonies without pigment or hemolysis at 48 hours on BAP; poor or no growth on MAC/EMB in 48 h; no distinctive odor

Reactions: Oxidase-variable; indole negative; catalase positive



Burkholderia pseudomallei (*Melioidosis*)

- Found in tropical areas
 - Southeast Asia & Australia
- Known as “**Vietnam time bomb**”
 - High levels in water of rice paddies–Vietnam & Thailand
 - Known to reactivate after many years of latency
- Types of infection:
 - Localized
 - Pulmonary
 - Bloodstream
 - Disseminated




Burkholderia pseudomallei

Collection and Transport of Clinical Specimens for Laboratory Rule-Out Testing	
Bone marrow or whole blood	<ul style="list-style-type: none"> • Considered the best specimen for culture • Collect directly into an appropriate blood culture bottle • Transport bottles at room temperature as soon as possible to obtain the diagnosis
Sputum or bronchoscopically obtained specimens	<ul style="list-style-type: none"> • Collect expectorated specimen into sterile transport cup or collect during bronchoscopy procedure. • Transport at room temperature up to 2h • If it is known that material will be transported from 2-24 h after collection, then store and transport at 2-8°C.
Tissue specimens (biopsies, abscess aspirates) and wound swabs	<ul style="list-style-type: none"> • Tissue pieces (at least the size of a pea) should be collected and kept moist • Transport in sterile container at room temperature within 1 hour of collection • Alternatively a swab from a tissue sample can be submitted in hospital transport tube with medium to stabilize specimen (e.g. Amies charcoal).
Urine	<ul style="list-style-type: none"> • Collect at least 1 mL into leak-proof container • Transport at room temperature up to 2 h • Refrigerate 2 up to 24 h until culture inoculation

Burkholderia pseudomallei

- Gram Stain
 - Straight, or slightly curved Gram negative rods
 - Colonies may demonstrate bipolar morphology in direct specimens and peripheral staining in older cultures, which can mimic endospores



Burkholderia pseudomallei

➤ Colony Morphology

- On BAP, small, smooth, creamy colonies in the first 1–2 days, that may gradually change in time to dry, wrinkled colonies (similar to *Pseudomonas stutzeri*)
- Poor growth at 24h, good growth at 48h
- Colonies are non-hemolytic and not pigmented on BAP. Organism grows on MAC
- Distinctive musty earthy odor
 - (the odor is apparent without sniffing or opening plate)


24 HOURS ON BAP



48 HOURS ON BAP



48 HOURS ON MAC



Burkholderia pseudomallei

➤ Biochemical Tests

- Oxidase positive
- Spot indole negative
- Motile
- Growth at 42 °C
- Polymyxin B and colistin no zone, penicillin resistant, amoxicillin–clavulanate susceptible

➤ Possible Misidentification:


- *Burkholderia cepacia*, *Chromobacterium violaceum*, *Pseudomonas aeruginosa*, *Pseudomonas stutzeri*, *S. maltophilia* and other nonfermenting Gram negative bacilli.

Major characteristics of *Burkholderia pseudomallei*:

Gram stain morphology: Gram negative rod, straight or slightly curved, may demonstrate bipolar morphology at 24 h and peripheral staining, like endospores, as cultures age

Colony morphology: Poor growth at 24 h, good growth of smooth, creamy colonies at 48 h on BAP; may develop wrinkled colonies in time, nonhemolytic. Can demonstrate strong characteristic musty, earthy odor; growth on MAC/EMB in 48 h, no pigment is visible on Mueller-Hinton agar, may have non-violet pigment on BAP.

Reactions: Oxidase positive; indole negative



Francisella tularensis (Tularemia)

- Highly virulent
 - Infectious dose 10–50
- Zoonotic–widespread throughout Northern Hemisphere
 - Hosts include: mammals, birds and amphibians
 - Primary host–small rodents
 - Has been reported in every state except Hawaii



Francisella tularensis

Forms of Tularemia

- **Ulceroglandular** This is the most common form of tularemia and usually occurs following a tick or deer fly bite or after handling of an infected animal.
- **Glandular** Similar to ulceroglandular tularemia but without an ulcer.
- **Oculoglandular** This form occurs when the bacteria enter through the eye.
- **Oropharyngeal** This form results from eating or drinking contaminated food or water.
- **Pneumonic** This is the most serious form of tularemia. Symptoms include cough, chest pain, and difficulty breathing.




Francisella tularensis

Collection and Transport of Clinical Specimens for Laboratory Rule-Out Testing	
<p>Tissue Biopsy, scraping of an ulcer, or conjunctival swab</p>	<ul style="list-style-type: none"> For small tissue samples, add several drops of sterile normal saline into a sterile container to keep the tissue moist Transport at room temperature for immediate processing If processing of specimen is delayed beyond 2 hours, keep specimen chilled (2-8°C) Amies transport media is an appropriate transport medium
<p>Aspirate Lymph node or lesion</p>	<ul style="list-style-type: none"> Submit in a sterile container Transport at room temperature for immediate processing. If processing of specimen is delayed beyond 2 hours, keep specimen chilled (2-8°C)
<p>Bone Marrow</p>	<ul style="list-style-type: none"> Submit in a sterile container Transport at room temperature for immediate processing If processing of specimen is delayed beyond 2 hours, keep specimen chilled (2-8°C)
<p>Blood</p>	<ul style="list-style-type: none"> Transport inoculated bottles directly to laboratory at room temperature Hold at room temperature until placed into automated, continuous monitoring blood culture incubators Do not refrigerate
<p>Respiratory Secretions</p>	<ul style="list-style-type: none"> Submit in a sterile container Transport at room temperature for immediate processing If processing of specimen is delayed beyond 2 hours, keep specimen chilled (2-8°C)
<p>Serum</p>	<ul style="list-style-type: none"> Collect at least 1 mL without anticoagulant for serologic diagnosis Store at 4°C until testing is performed. Acute specimen is collected as soon as possible after onset of disease Convalescent-phase should be collected >14 days after the acute specimen.

Francisella tularensis

- Gram Stain:
 - Tiny, Gram negative coccobacilli
 - Poorly counterstaining with safranin (basic fuchsin counterstain may increase resolution)
 - Pleomorphic
 - Mostly single cells




Francisella tularensis

➤ Colony Morphology

- Aerobic, fastidious
- No growth on MAC or EMB
- Scant or no growth on BAP; may grow on primary culture, not well on subculture
- Slow growing on CHOC, TM or BCYE: 1–2 mm after 48h
- Colonies are opaque, grey–white, butyrous with smooth and shiny surface

48 hours on BAP


- After 48 h: <1 mm, gray-white, opaque, no hemolysis.



Francisella tularensis

➤ Biochemical Tests

- Oxidase negative
- Catalase negative or weakly positive
- Satellite negative
- Beta-lactamase positive




F. tularensis is beta-lactamase positive (+).

- Common Misidentifications: *Acinetobacter*, *Aggregatibacter actinomycetemcomitans*, *Haemophilus influenzae*, *Oligella* spp. and *Psychrobacter* spp

Major characteristics of *Francisella tularensis*

Gram stain morphology: Pleomorphic, minute (0.2 to 0.5 by 0.7 to 1.0 μm) faintly staining, Gram negative coccobacillus


Colony morphology: No growth on MAC/EMB, scant to no growth on BAP after > 48 h. Produces 1-2 mm gray to grayish-white colonies on CHOC after > 48 h



Yersinia pestis (Plague)

- Plague is a zoonotic disease
- There are 1–17 human cases per year in the U.S. (50% are from New Mexico)

Reported cases of human plague--United States, 1970-2012



1 dot placed in county of exposure for each plague case

Yersinia pestis

➤ Host range:

- A number of species of rodent fleas act as natural vectors.
- Rodent reservoir (squirrels, prairie dogs, chipmunks, deer mice, voles, rats, occasionally hares & rabbits).
- Most mammals, including humans and carnivores (such as bobcats, mountain lions, dogs & cats) are accidental hosts; not birds or reptiles.
- Cats & Dogs





Yersinia pestis

➤ **Bubonic**

- Transmitted by the bite of an infected flea
- Untreated mortality ~60%,
- Not transmitted person-to-person

➤ **Septicemic**

- Results from contact with, or being bitten by, an infected animal
- Petechiae, necrosis and gangrene of the extremities (“the black death”)
- Untreated mortality ~100%
- Not transmitted person-to-person

➤ **Pneumonic**

- Inhalation– deadliest form–intentional release as an aerosol
- May be transmitted person-to-person, or animal-to-person, extremely contagious

Yersinia pestis

Bubonic Plague



Pneumonic Plague



Septicemic Plague




Yersinia pestis

Collection and Transport of Clinical Specimens for Laboratory Rule-Out Testing	
Lower respiratory tract	<ul style="list-style-type: none"> • Transport specimens in sterile, screw-capped containers at room temperature. • If it is known that material will be transported from 2-24 h after collection, then store and transport at 2-8°C
Blood	<ul style="list-style-type: none"> • Transport samples directly to the laboratory at ambient temperature and place onto the blood culture instrument • Do not refrigerate • Follow established laboratory protocols for processing blood cultures
Aspirate, tissue or biopsy specimen	<ul style="list-style-type: none"> • Submit tissue or aspirate in a sterile container. • For small samples, add 1–2 drops of sterile normal saline to keep the tissue moist. • Transport sample at room temperature for immediate processing. • Keep the specimen chilled if processing will be delayed (> 2 h).
Swabs	<ul style="list-style-type: none"> • A swab of tissue is not recommended. • However, if a swab specimen is collected, the swab should be reinserted into an appropriate transport package and sent to the laboratory at room temperature for immediate processing. Keep the specimen chilled if processing will be delayed (>2 h).

Yersinia pestis


➤ Gram Stain

- Plump Gram negative rods seen mostly as single cells or pairs, and may demonstrate short chains in liquid media
- May exhibit bipolar, “safety-pin” appearance that is not seen on Gram stain, may be exhibited by Giemsa stain or Wrights stain




Gram stain of *Yersinia pestis* from a patient blood culture bottle

Note chaining and bipolar-staining Gram negative rods



Yersinia pestis

- Colony Morphology
- Facultative anaerobe
- Slow growing at 35 °C, better growth at 25–28 °C
- Grey–white, translucent pinpoint colonies at 24h
- At 48 h, colonies on BAP gray–white to slightly yellow and opaque
- “Fried egg” on BAP in older cultures (96h)
- At 48h, lactose non–fermenter on MAC or EMB



48 HOURS ON BAP at 25°C and 35°C

Yersinia pestis

➤ Biochemical Tests

- Catalase positive
- Oxidase, urease (at 35 ° C) and indole: negative

➤ Common Misidentifications:

- *Shigella* spp., H₂S(-) *Salmonella* spp., *Acinetobacter* or *Pseudomonas* spp., or *Yersinia pseudotuberculosis*.
- May not be identified in common automated ID systems, including MALDI TOF

Major characteristics of *Yersinia pestis*:

Gram stain morphology: Gram negative rods, 0.5 x 1-2 μm

Colony morphology: Slow growing at 35°C with either pinpoint colonies or no growth on BAP after 24h; colonies are 1-2 mm, gray-white to slightly yellow and opaque on BAP after 48 h; non-lactose fermenter on MAC/EMB

Specimen is blood, sputum, or lymph node aspirate

Yes

Oxidase: Negative
Catalase: Positive
Indole: Negative
Urease: Negative

No

Not *Yersinia pestis*.
Continue identification per
routine laboratory procedures.
May be other *Yersinia* spp.

Yes

Y. pestis is not ruled out.

Contact LRN Reference Level Laboratory to refer the isolate.

Report: *Y. pestis* cannot be ruled out; isolate referred to LRN Reference Level Laboratory.



Next Steps for Select Agent Identification




Select Agent ID

- Next steps once you cannot rule out a select agent:
 1. Notify the LRN (MDHHS Bureau of Laboratories)
 2. Properly package and ship isolate to MDHHS for further testing– (i.e. Category A)
 3. Secure all related specimens and isolates until confirmatory testing is complete by the LRN

Select Agent ID

SELECT AGENT ALGORITHM GUIDE






Select Agent ID

- Next steps once MDHHS confirms select agent
 1. MDHHS will contact sentinel lab with results
 2. MDHHS Lab/EPI will discuss potential lab exposures
 1. Exposed lab staff may require monitoring & prophylaxis
 2. If exposures occurred–CDC form 3 will need to be completed
 3. MDHHS will notify the CDC–CDC form 4 is completed
 1. MDHHS fills out section A & B of CDC form 4
 2. Sentinel Lab fills out section C & D of CDC form 4
 4. Sentinel lab destroys all isolates, specimens and derivatives within 7 days of MDHHS select agent confirmation

Select Agent ID

SELECT AGENT ALGORITHM GUIDE





Automated Systems Study



Preparation Safety Study

- Determined viability of BT agents using direct, extended direct, and tube extraction
 - *Bacillus anthracis* Sterne, *Brucella abortus* Strain 19, *Francisella tularensis* LVS, *Yersinia pestis* A1122, *Burkholderia thailandensis* ATCC 70038, *Clostridium* spp (*botulinum* types A, B, E, and *perfringens*)
- Some isolates of *B. anthracis*, *C. botulinum*, *F. tularensis*, *Y. pestis*, and *B. abortus* survived the direct and extended direct prep methods
- None survived the tube extraction
- *B. anthracis*, *B. thuringiensis*, *B. cereus*, *B. subtilis* survived tube extraction in other studies
- Use tube extraction AND 0.2 μm filter for “hazardous organisms”

Automated Systems Study

- Determine ability of Bruker and Vitek software libraries to accurately identify BT agents
 - Test IVD, ROU, and Security Libraries independently
- Panel of 6 strains and 6 near neighbors for each agent
 - *B. anthracis*, *Y. pestis*, *F. tularensis* (5), *C. botulinum* (4), *B. mallei*, *B. pseudomallei*, *Brucella* spp.
 - Select isolates from AOAC Stakeholder Panel on Agent Detection Assays (SPADA) – inclusivity and exclusivity panels
 - Each isolate tested in triplicate
 - Acceptable ID:
 - Bruker sample score ≥ 2.0 ; Vitek ID score $\geq 60\%$
- Challenges
 - No SPADA panels for *C. botulinum* or *Brucella* spp.
 - Isolate availability and shipping
 - No tube extraction procedure for Vitek MS



Study Participants

- Michigan Department of Health and Human Services
- New York Department of Health, Wadsworth Center, Biodefense Laboratory
- North Carolina State Laboratory of Public Health
- New York City Department of Health and Mental Hygiene
- Florida Department of Health, Bureau of Public Health Laboratories – Jacksonville Branch
- Minnesota Department of Health
- State Hygienic Laboratory at the University of Iowa
- Texas Department of State Health Services, Laboratory Services Section, Austin
- Center for Microbial Genetics and Genomics, Northern Arizona University

Study Limitations

- Samples were prepared by four labs, shipped and stored frozen
- Some participants tested extracts beyond the 45 day time frame established in the study protocol
- Some participants prepared new targets for testing each software library

Organism (# of strains tested)	Reported Identification	IVD Library		RUO Library		Security Library	
		Mean Score (# replicates)	Species level ID % (>2.0)	Mean Score (# replicates)	Species level ID % (>2.0)	Mean Score (# replicates)	Species level ID % (>2.0)
<i>B. anthracis</i> (6)	No reliable ID	1.13 (90)		1.38 (45)		1.1 (40)	
	<i>B. cereus</i>			1.88 (62)	8.3		
	<i>B. anthracis</i>					2.08 (68)	50
	<i>B. pseudomycooides</i>			1.56 (1)			
<i>B. thuringiensis</i> (1)	No reliable ID	1.1 (12)		1.33 (9)		1.12 (9)	
	<i>B. cereus</i>			2.01 (9)	22.2		
	<i>B. anthracis</i>					2.06 (9)	38.9
<i>B. circulans</i> (1)	No reliable ID	1.06 (12)		1.33 (5)		1.01 (18)	
	<i>B. circulans</i>			1.87 (13)	16.7		
<i>B. cereus</i> (1)	No reliable ID	1.15 (12)		1.35 (9)		1.23 (9)	
	<i>B. cereus</i>			2.08 (9)	44.4		
	<i>B. anthracis</i>					2.14 (9)	50
<i>B. mycooides</i> (1)	No reliable ID	1.15 (15)		1.61 (7)			
	<i>B. mycooides</i>			1.77 (7)			
	<i>B. anthracis</i>					1.77 (18)	
	<i>B. weihenstephanensis</i>			1.19 (4)			
<i>B. megaterium</i> (1)	No reliable ID	1.12 (12)		1.35 (3)		1.06 (18)	
	<i>B. megaterium</i>			1.98 (15)	55.6		
<i>B. subtilis</i> (1)	No reliable ID	1.04 (15)		1.47 (13)		0.98 (18)	
	<i>B. subtilis</i>			1.68 (5)			

Organism (# of isolates tested)	Reported Identification	IVD Library		RUO Library	
		Mean Score (# replicates)	Species level ID % ($\geq 60\%$)	Mean Score (# replicates)	Species level ID % ($\geq 60\%$)
<i>Y. pestis</i> (6)	No Identification	0 (2)		0 (23)	
	<i>Y. pseudotuberculosis</i> / <i>Y. frederiksenii</i>	51.7/48.2 (7)			
	<i>Y. ruckerii</i> / <i>Y. pseudotuberculosis</i> / <i>Y. frederiksenii</i>	33.3/33.4/33.2 (1)			
	<i>Y. pseudotuberculosis</i>	99.9 (17)	58.6	87.8 (18)	33.3
	<i>Y. pestis</i>	99.9 (1)	3.4		
	<i>Yersinia</i> spp.			88.8 (13)	
<i>Y. ruckeri</i> (1)	No Identification	0 (5)			
	<i>Y. ruckerii</i>	99.7 (4)	44.4		
	<i>Yersinia</i> spp.			85.7 (9)	
<i>Y. pseudotuberculosis</i> (3)	No ID	0 (1)			
	<i>Y. pseudotuberculosis</i>	99.6 (20)	95.2	91.9 (20)	95.2
	<i>Y. enterocolitica</i>			94.1 (1)	4.8
<i>Y. enterocolitica</i> (2)	No Identification	0 (1)			
	<i>Y. enterocolitica</i>	99.9 (8)	66.7	97.5 (11)	91.7
	<i>Y. pseudotuberculosis</i>			86 (1)	
	<i>Y. pseudotuberculosis</i> / <i>Y. enterocolitica</i>	49.1/50.9 (3)			



Identification Errors

Organism	Concern	Organism	Concern
<i>B. anthracis</i>	<i>B. cereus</i> <i>B. cereus</i> group Split identifications	<i>B. mallei</i>	<i>B. pseudomallei</i>
<i>B. thuringiensis</i>	<i>B. anthracis</i>	<i>B. pseudomallei</i>	<i>B. mallei</i>
<i>B. cereus</i>	<i>B. anthracis</i>	<i>B. thailandensis</i>	<i>B. pseudomallei</i> <i>B. mallei</i>
<i>Y. pestis</i>	<i>Y. pseudotuberculosis</i> <i>Yersinia</i> spp. Split identifications	<i>C. botulinum</i>	<i>C. sporogenes</i>
<i>Y. pseudotuberculosis</i>	<i>Y. pestis</i>	<i>B. melitensis</i>	<i>Brucella</i> spp.
<i>Y. enterocolitica</i>	<i>Y. pestis</i>	<i>B. abortus, suis, canis, ovis, ceti pinnipedialis</i>	<i>B. melitensis</i> <i>Brucella</i> spp.
<i>F. philomiragia</i>	<i>F. tularensis</i>		



Study Summary

Software libraries for BT agents need improvement

- **Currently not reliable for BT agents**
- Addition of BT agents to IVD software
- Availability of alternative libraries i.e. MicrobeNet

Laboratories must be aware of software limitations

- Partial or incorrect identifications
- Additional training by manufacturer or BT coordinator

Continue use of the rule out and refer protocols

