

Safety made simple.

BD Eclipse™ Needles. Important Usage Guidelines.

Step 1.



Firmly attach needle onto a Luer Lock* syringe with a push and clockwise twist.

Step 2.



Pull back on the safety cover toward the syringe and away from the needle.

Remove needle shield.

Step 3.



Draw up medication as per your usual technique.

NOTE: If changing needles, activate safety cover prior to removing the needle from the syringe.



Inject according to protocol.

For user convenience, the needle "bevel-up" position is oriented to the safety cover.

*For use with Luer Lock syringe only.

Step 4.



After injection, immediately activate safety cover as follows:

- Using the same hand holding the syringe, center thumb or forefinger on the finger pad area of the safety cover.
- Push the cover forward toward the needle until you hear and/or feel it lock.



Visually confirm that the needle tip is covered.

After activation, discard immediately into an approved sharps container.

NEW ISO NEEDLE HUB/PACKAGING CHART

21G	25G	23G	22G	18G
-----	-----	-----	-----	-----

BD Medical Systems
1 Becton Drive
Franklin Lakes, NJ 07417

www.bd.com