

# Needlestick prevention that meets your needs



- virtually eliminates exposure
- single-handed, passive activation
- hands remain behind the sharp
- integrated safety mechanism
- prevents recapping and needle removal
- remains safe through disposal
- easy to use, requires minimal training
- safe and effective for patients

Catalog Number	Description
<b>VanishPoint® Allergy Syringe Tray</b> 25 per tray/1000 per case	
10134	1ml VanishPoint allergy syringes 27G x 1/2"
<b>1ml VanishPoint® syringes</b> 100 per box/800 per case	
10131	1ml VanishPoint tuberculin syringe 27G x 1/2"
10151	1ml VanishPoint tuberculin syringe 25G x 5/8"
10161	1ml VanishPoint tuberculin syringe 25G x 1"
10211	1ml VanishPoint U-100 insulin syringe 29G x 1/2"
<b>3ml VanishPoint® syringes</b> 100 per box/600 per case	
10301	3ml VanishPoint syringe 25G x 5/8"
10391	3ml VanishPoint syringe 25G x 1"
10311	3ml VanishPoint syringe 23G x 1"
10321	3ml VanishPoint syringe 23G x 1 1/2"
10331	3ml VanishPoint syringe 22G x 1"
10341	3ml VanishPoint syringe 22G x 1 1/2"
10351	3ml VanishPoint syringe 21G x 1"
10361	3ml VanishPoint syringe 21G x 1 1/2"
10371	3ml VanishPoint syringe 20G x 1"
10381	3ml VanishPoint syringe 20G x 1 1/2"
<b>5ml VanishPoint® syringes</b> 100 per box/600 per case	
10531	5ml VanishPoint syringe 22G x 1"
10541	5ml VanishPoint syringe 22G x 1 1/2"
10551	5ml VanishPoint syringe 21G x 1"
10561	5ml VanishPoint syringe 21G x 1 1/2"
10571	5ml VanishPoint syringe 20G x 1"
10581	5ml VanishPoint syringe 20G x 1 1/2"
<b>10ml VanishPoint® syringes</b> 100 per box/600 per case	
11031	10ml VanishPoint syringe 22G x 1"
11041	10ml VanishPoint syringe 22G x 1 1/2"
11051	10ml VanishPoint syringe 21G x 1"
11061	10ml VanishPoint syringe 21G x 1 1/2"
11071	10ml VanishPoint syringe 20G x 1"
11081	10ml VanishPoint syringe 20G x 1 1/2"
<b>VanishPoint® blood collection products</b>	
22701	VanishPoint blood collection tube holder 250 units per case
22711	VanishPoint small diameter tube adapter 25 units per box



**RETRACTABLE  
TECHNOLOGIES, INC.**  
www.vanishpoint.com

© 2006 Retractable Technologies, Inc.  
All rights reserved.  
M005 REV. 5

511 Lobo Lane  
Little Elm, Texas 75068-0009 USA  
Tel: (972) 221-6644  
Fax: (972) 294-4400  
Toll Free: 1-888-703-1010

domestic: rtisales@vanishpoint.com  
international: rti.intl@vanishpoint.com

# UNCOMPROMISING SAFETY FOR CLINICIANS AND PATIENTS



**VANISHpoint™**  
SYRINGES

The New Standard for Safety™

# VANISHPOINT® Safety by design

- Safety mechanism is built into the syringe
- Needle retracts directly from patient into syringe, effectively reducing exposure to contaminated sharps
- Attached needle design:
- maintains a closed system
  - prevents leakage and contamination
  - reduces deadspace
  - prevents unnecessary needle changes
- "Changing needles between drawing vaccine from a vial and injecting it into a recipient is unnecessary."\*
- Safe for you and your patients



\*CDC (2002). General recommendations on immunization: recommendations of the Advisory Committee on Immunization Practices (ACIP) and the American Academy of Family Physicians (AAFP). MMWR 51(RR02):1-36.

1

Use standard procedure for drawing up medication.



2

Prepare and give injection using aseptic technique according to institutional policy.



3

While the needle is still in the patient, fully depress the plunger to activate retraction.



4

Dispose of VanishPoint® syringe in an appropriate sharps container.



# the **VANISHPOINT®**

## **ADVANTAGE**

**Designed for Clinician and Patient Safety**

Clear, unobstructed calibrations allow for accurate dosing



Triple beveled, lubricated needle provides patient comfort



**NO CONTAMINATED SHARP EXPOSED**

Attached needle prevents leakage and contamination

Pre-removal activation prevents exposure to contaminated sharp



**REDUCED DEADSPACE**

VanishPoint® syringes are available with a variety of needle sizes and gauges attached. Syringes are color coded for needle gauge.

VanishPoint® tube holders are used with standard blood collection needles and vacuum tubes for safe blood collection



*Activated VanishPoint® Syringes require less disposal space than other syringes and prevent disposal-related injuries.*



VanishPoint® Allergy Syringe Tray