September 11, 2014



# Update #1 - Enterovirus D68 (EV-D68)

**Target audience:** Clinical laboratories

## **Background:**

See Enterovirus D68 (EV-D68) Guidance 9-8-2014
 http://www.michigan.gov/documents/mdch/Enterovirus+D68+Guidance+9-8-14 467990 7.pdf

# **Current Status:**

- CDC has not issued formal guidance about ancillary testing for suspect patients. However, due to a high number of test requests, it is suggested that testing for other agents associated with severe respiratory illness be conducted.
- Testing should include: respiratory syncytial virus (RSV), influenza, and rhino/enterovirus.

## **Laboratory Testing for Enterovirus D68:**

Authorization for testing is required. Please contact MDCH Bureau of Disease Control, Prevention and Epidemiology at: 517-335-8165. After Hours call: 517-335-9030

- Currently, testing is being performed by CDC.
- Specimens should be submitted to MDCH Bureau of Laboratories and will be shipped to CDC.

# **Preferred Specimens for Enterovirus D68:**

- NP swabs (preferred) in viral transport
- ½ 1ml Serum
- Cerebral spinal fluid, if neurological symptoms are present
- Store at 4°C

### Packaging and Shipping Clinical Specimens to MDCH BOL

**Please complete both the MDCH Specimen and CDC Submission Forms.**All specimens should be packaged and shipped to MDCH BOL on cold packs

Link to MDCH Specimen Submission Form:

http://www.michigan.gov/documents/DCH-0583TEST\_REQUEST\_7587\_7.pdf

(Write or type the test requested, i.e. "Rule-out EV-D68" in blank space found under "Hepatitis" at the bottom right of page 1.)

Link to CDC Specimen Submission Form:

http://www.michigan.gov/documents/mdch/HUMAN form-50-34 410210 7.pdf

Tips for completing the CDC form:

• Test Code: CDC-10312

• Test Order Name: Enterovirus detection and Identification

Suspect agent: LEAVE BLANK

Comments on page 2: Rule-out EV-D68

### Steps for laboratories to take now:

- Please monitor the situation and be aware of current guidance.
- Share this information with others in your facility.

#### Questions and additional information:

For laboratory related questions, please contact:

Dr. Janice Matthews-Greer, Virology Section Manager

<u>MatthewsGreerJ@michigan.gov</u> (517) 335-8063 (during normal business hours)