

MICHIGAN DEPARTMENT OF HEALTH AND HUMAN SERVICES

BUREAU OF LABORATORIES

Electronic Test Ordering and Results (ETOR) Training & User Manual

TABLE OF CONTENTS

Introduction	3
What is ETOR? What isn't ETOR?	3
The LifeCycle of a Sample in ETOR	3
SETTING UP YOUR STATE OF MICHIGAN SINGLE SIGN ON ACCOUNT	4
REQUESTING ACCESS TO ETOR	13
SIGN UP FOR ETOR APPLICATION FOR SPECIMEN DATA ENTRY	18
LOGGING INTO ETOR	23
Data Entry for New Samples	26
Packaging and Submitting Samples	29
Multiple Clinic Mary	32
NAVIGATION BAR	32
SEARCH/PACKAGE SCREEN	33
Editing Test Orders (Before Packaged)	34
Canceling Test Orders (Before Packaged)	34
Packaging and Shipping Your Samples	35
LOOKING UP RESULTS	37
AN EXAMPLE OF A FINAL REPORT	38
AVAILABLE TESTS	39
THE FOLLOWING TESTS ARE CURRENTLY AVAILABLE IN ETOR:	39
WHERE TO GET HELP	40

Introduction

What is ETOR? What isn't ETOR?

ETOR stands for Electronic Test Ordering and Results system.

ETOR is

 A portal to the State Laboratory Information System for Submitters to enter demographics on specimens being submitted to the MDHHS Bureau of Laboratories, and view/download results for samples that have completed testing in the MDHHS Laboratory and were submitted via ETOR.

The ETOR portal provides:

- A view of sample information
- Package slips to expedite receipt at MDHHS Laboratory
- The ability to track the sample as it flows through the MDHHS Laboratory
- The ability to view and download test results (via pdf file) after testing is completed.

ETOR isn't

- Involved in HL7 messaging
- Designed for the Health Interoperability Exchange programs

The LifeCycle of a Sample in ETOR

A sample submitted via ETOR has a lifecycle as follows.

Sample Status:


- 1. (NONE) Pre-sample data entry. Once it is saved, it will be "IN PROCESS."
- 2. **IN PROCESS** Sample(s) are entered into the ETOR portal.
 - a. A sample's submission data is saved in ETOR. It will appear on the ETOR Search/Package Screen.
 - b. A sample(s) data can be edited within ETOR until they are packaged.
 - c. In Process orders can be canceled. The sample status will change to "CANCELED".
- 3. **SUBMITTED** Sample(s) have been selected on the ETOR Search/Package Screen, and have been Packaged.
 - a. A package slip will be available for printing to ship to the laboratory with the sample(s). A package slip can contain one or more specimens.
 - b. At this point a sample's data can no longer be edited in ETOR and the order can no longer be canceled
- 4. **SUBMITTED ACK** The Laboratory's Information Management System (LIMS) has electronically acknowledged and received the information provided for the sample(s) packaged in ETOR.
 - a. The package of specimen(s) has not yet been received by the Laboratory, but the LIMS has a copy of the sample(s) information submitted in ETOR.
- 5. **RECEIVED** The data entry staff at the MDHHS Laboratory have received the package and scanned the barcode on the packing slip. The sample is in the LIMS and is ready to be tested.
- 6. **COMPLETED** Testing at the MDHHS Laboratory has been completed and the final result is viewable and downloadable via the ETOR Search/Package Screen.

Setting Up Your State of Michigan Single Sign On Account


If your facility already has a Single Sign On account, use your current log-in and skip to "Requesting Access to ETOR" below.

If your facility does not already have a Single Sign On account, navigate to the following URL-https://sso.state.mi.us/som/dch-portal/dch-portal


- 1. Click the "Register" button.
- 2. Fill out the form. Be certain your email address is correct.




3. Select a four digit number to create your unique UserID and enter the number in the box on the form.



4. Review your information. If incorrect, go back to edit. Otherwise click "Submit."



5. You have now applied for a Single Sign On account. After approval, you will need to request access to the ETOR application link.




6. Once approved, follow the "All Other Application Users click here to login" link.




If you are a new Single Sign On account user and have any problems accessing your account, please contact the State of Michigan's Client Service Center at 241-9700 or 1-800-968-2644.

Michigan Business One Stop clients only - If you have any problems accessing your One Stop account please contact the Customer Assistance Center at 877-766-1779.


7. Login with your username and temporary password.




8. Change your password using the Single Sign On password rules.



9. Enter your answers to the security questions.




10. Click "Done."




Requesting Access to ETOR


- 1. Log into the Michigan Single Sign On application.
- 2. To subscribe to the MDHHS Laboratory ETOR application, click on the "Subscribe to Applications" link.




3. Select the "Dept of Community Health" from the drop down list and "MDCH Laboratory Electronic Test Ordering and Results" from the "Select App" drop down.




4. Enter a work phone number to be used if the MDHHS Laboratory needs to confirm your approval for ETOR. The email that you used to create a SSO account will be automatically populated into the email box.



5. Confirm your information.




6. You are now subscribed and will be able to see the ETOR link on the Single Sign On page. Next you will need to request approval in the ETOR portal to begin logging in samples.






Sign Up for ETOR Application For Specimen Data Entry

- 1. Sign on to your Single Sign On account.
- 2. Click on the "MDCH Laboratory Electronic Test Ordering and Results" link to begin.




3. Accept the MDCH Systems Use agreement.




- 4. Enter your Submitter Agency information. If you don't know your StarLIMS agency number you can:
 - a. Find your agency number on the final reports you receive from the MDHHS Bureau of Laboratories. This can be located in the report header along with the submitter address.
 - b. Contact the Data and Specimen Handling Unit via fax (517-335-9871) on your agency letterhead and request your agency's "Submitter ID."

NOTE: Entering your StarLIMS agency number will help expedite your user approval. NOTE: The same submitter may be associated with more than one StarLIMS agency number. This is likely when your agency submits samples for lead AND clinical (ex. CTGC, HIV etc) testing (note that the StarLIMS agency number for lead submitters is alpha-numeric ex.AAA##). Please see the

Muliple Clinic Mary section of this manual if this applies to your agency.



5. Verify your information is correct.




6. Once your application is approved, you will get a confirmatory email like the one below. It may take a few business days for the MDHHS Laboratory to approve your ETOR use. Your application may be rejected if the MDHHS Laboratory cannot confirm which facility you are requesting access to or if your facility does not verify your employment and need for access.




Logging into ETOR

1. Log into your State of Michigan Single Sign On account and click on the "MDCH Laboratory Electronic Test Ordering and Results" link.



2. Accept the MDCH Systems Use Agreement.



3. ETOR Main Menu.

There are two options on this screen. One enables tests to be ordered on new specimens and the other allows a search for previously entered orders or allows previously entered specimens to be packaged for shipping.




Data Entry for New Samples


- 1. Click the "Add New Test Order' button.
- 2. Select the Laboratory you intend to send the specimen to and the test(s) you are requesting for that sample. You may select multiple tests, but only one Laboratory.

NOTE: You should only order multiple tests against an order if the same **physical** sample is being used to complete those tests (1 ETOR Order ID = 1 sample). For example, if Patient A requires Syphilis (USR) and CTGC testing and you have collected a serum and urine sample then two separate orders need to be created in ETOR. On the other hand, if Patient A requires Syphilis (USR) and HIV Ag/Ab Serum and you have collected a serum sample then both of these tests can be ordered against the same ETOR Order ID.


3. Click "Continue" to view the Test Order / Data entry screen.



- 4. Fill in the on-line form as completely as possible.
- 5. Required fields are marked with a red asterisk, and the sample cannot be saved until all required fields are filled in. The "Specimen ID" field is used for as an ETOR searchable identifier to quickly locate this sample in the system. Anything placed in the "Comments" field are visible to the MDHHS Laboratory Data Entry Staff. If you would like the results to be cc'd to another agency, a comment can be entered to indicate this (note that, as an alternative, this comment could also be hand-written on the package slip).




6. Each sample will receive a unique "Order ID."



Packaging and Submitting Samples

- 1. In the view/edit package screen all orders for your agency can be located. The "Status" field indicates the status of each sample: None, In Process, Submitted, Submitted ACK, Received, Completed, and Canceled. (For a review of what each status indicates, see the section titled <u>The LifeCycle of a Sample in ETOR</u> of this manual.)
- 2. In this example, the sample hasn't yet been packaged as indicated by the STATUS indicator "In Process."




- 3. To package it, simply check the box to the left of all samples to be packaged and click "Create Package." This will create a package identifier, and pop up a printable form.
- 4. Example of the printable package slip:



NOTE: There is a barcode associated with the Package ID and each Order ID as shown above. If you notice a number instead of the barcode (s) when displayed/printed, please try using another browser. Internet Explorer is recommended.

5. Going back to the view/edit package screen, note that the specimen status has changed to "Submitted," indicating it's ready to mail to the MDHHS Laboratory. Please remember to print and include the above packing slip in each box of samples being submitted to the Laboratory for testing. For more information on packaging and submitting samples, see the **Search/Package Screen** section of this manual.



Multiple Clinic Mary

"Multiple Clinic Mary" describes the situation of a nurse who works at two clinics. At each clinic the nurse, Mary, is not allowed to see confidential HIPAA data from the other clinic. The web application tracks which sites a user is logging in from. When "Mary" logs on at one site, the system first uses IP addresses to determine which site Mary is logging on from. It then checks to ensure that Mary is authorized to work at that site. Once Mary logs in, she can see only patient data from her current site. In this situation, Mary will **need to have multiple logins/passwords, one for each site**.

Navigation Bar

The navigation bar is present on all screens *except* the login page. The navigation bar gives basic links, such as:

- The Michigan.gov home page
- ETOR home screen / main menu
- Contacting ETOR Administrators
- MDCH Home Screen
- ETOR Manual
- Log out



Search/Package Screen

From this page, a user can see all the tests that were entered, their status, whether they were packaged / submitted to the lab, and the final result if testing has been completed by the laboratory.

From this screen you can:

- a) Edit or cancel the test order **before** it has been packaged
- b) Records Per Page: Click the down arrow (highlight number of records per page)
- c) Cycle through pages of the view
- d) Search orders
- e) See order status (See "The Life Cycle of a Sample in ETOR" section)
- f) Sort selections: Click once to descend, click again to ascend the sort selection fields. (Date Created, Collection Date, Patient ID, etc.)
- g) Select order(s) to be packaged generating a package slip
- h) View a package slip



NOTE: If the Search/Package screen grid extends off the right side of the page (as shown below) the following can be done to correct this: If you are in Internet Explorer and you click on Tools, your "Compatibility View" should be unchecked. If it is already unchecked, then go to Tools, Compatibility View Settings, and make sure that the "Display intranet sites in Compatibility View" and "Display all websites in Compatibility View" are unchecked.




Editing Test Orders (Before Packaged)

To edit information for an order, click on the view/edit order # in the Search/Package screen. This will return you to the Test Order screen for this sample. If it has not yet been packaged, you will be able to edit and save the information. See section below for instructions on how to deal with errors after an order was packaged.

Canceling Test Orders (Before Packaged)

Orders can be canceled if they are in a status of IN PROCESS. To cancel an order(s), check the box to the left of all orders to be canceled in the Search/Package screen and click "Cancel Orders(s)".



If successful, the notification "Order(s) successfully canceled" will appear and the status will be updated to CANCELED.

NOTE: If an order is canceled by mistake you can select the order and click "Create Package" to submit the order and continue with the usual workflow.

Packaging and Shipping Your Samples

Samples submitted to the MDHHS Laboratory must comply with federal regulations to ensure testing accuracy. Therefore, all samples submitted to MDHHS must adhere to federal packaging and labeling standards.

Be sure to observe these points when submitting your samples:

- For CT/GC, ship no more than 5-7 samples together in one package. For Rabies, each specimen must have its own package slip (1 order/slip) to allow the specimens to be distinguished. For other tests, you may ship any number of samples per package. All requests must include the ETOR Packing Slip with its corresponding samples. ETOR calls each set of samples with an accompanying cover sheet a "package." You can include multiple packages together in a single shipment, but each package should be wrapped separately, with its corresponding Packing Slip sheet.
- Each tube must be labeled with the patient's name, unless it is an anonymous HIV tube. The name on the tube must exactly match the name on the ETOR coversheet. **Handwritten corrections to the cover sheet are not acceptable.** For anonymous HIV samples with no names, a unique identifier should be used rather than patient name. In this case, the unique identifier on the tube must exactly match the unique identifier on the cover sheet. **Samples where the tube and cover sheet do not match will be rejected.**
- Pre-paid CT/GC samples can be ordered through ETOR by including the pre-paid form in the shipping box with the samples. The only information needed to be hand written on these pre-paid forms is the patient name.
- If you find an error in the package before you click "Create Package" you can edit the tests and fix the incorrect information or cancel the order. If you find an error after you package the orders, you should re-enter the incorrect test(s), print a new package slip and send us both package slips (the one with the error and the one without). Hand write a comment on the erroneous package slip so that the lab sample receiving unit is aware of the error. The erroneous order will be resulted as Test Not Done so that it is cleared from the system. The same procedure should be followed for cases where duplicate orders are created in ETOR by mistake.
- If you would like the results to be cc'd to another agency, a comment can be entered in the "Comments" field on the Data Entry Screen to indicate this. As an alternative, this comment could also be hand-written on the package slip.

To Print a Package Slip

- A. Click: View Package Slip #### (to view the package slip of choice)
- B. Right Click (anywhere on the package slip) and click Print.

When the user clicks **Create Package**, a package is created the selected orders are marked as belonging to that package. Those orders are now eligible for electronic acknowledgement by the laboratory system and their status will change to reflect this.


At this point, a package slip will pop up displaying information about the created package. The package slip contains information including barcodes that the Laboratory will use when receiving the package.

Note that browsers with popup blocking software will prevent the package slip from appearing. If this occurs, you can also view the package slip by going to the Package Status screen (see below).

To prevent your browser from blocking the package slip from popping up, you can do one of the following options, depending on what browser and popup blocking software you are running:

- Disable your popup browser software
- Hold the Control key down when submitting the package
- In Internet Explorer, select Tools from the menu at the top of the screen, then choose Popup Blocker, then Popup Blocker Settings. Enter this site in the field "Address of Web site to allow."

At this time the Firefox browser is not capable of printing the barcodes. Please use an alternate browser to enable package slip printing.



Looking Up Results

The Search/Package screen allows the user to find out whether the test has been completed and view or download the preliminary or final result report.


To look up a test submission, click the Search button on the lower left side of the screen. Select field(s) to search by. Enter data to search on – order number, package ID number, DOB, Specimen ID, etc. Once you have located your order, click the View Report link, and print or save the result PDF.

To view a result report, select the View Results link for that order. If there is no View Results link available, a report has not been associated with this sample. If the Status is COMPLETE or PRELIMINARY but no View Results link is available, contact the lims_help@michigan.gov team via email. Be sure to have the Subject line of your email start with "ETOR." Include the order number(s) affected and the date and time you are seeing the issue.

The number of reports displayed per page can be changed by using the "Records per Page option at the bottom. To navigate through additional reports per page select the numeric page menu at the bottom of the report list.

An Example of a Final Report

Below is a sample report.



Available Tests

The following tests are currently available in ETOR:

- Aerobic Isolate ID Environmental
- AFB Identification Isolate ID
- AFB Nucleic Acid Amplification *
- AFB Slide/Culture Clinical Specimen
- AFB Slide/Culture Clinical Specimen Non-Human
- Arbovirus Encep. Panel (IgM) *
- Bacterial Respiratory PCR Panel
- C. trachomatis & N. gonorrhoeae Non-Culture
- C. trachomatis Non-Culture
- E. coli Cultural Isolate Environmental
- E. coli (SLT) Toxin & Serology
- EIA STX
- Enteric Bacterial Culture
- Enterovirus PCR *
- Fungal Identification Isolate ID
- Fungal Serology Complement Fixation
- Hepatitis A Antibody (IgM)
- Hepatitis B Antibody (Anti-HBsAg)
- Hepatitis B Surface Antigen (HBsAg)
- Hepatitis B Surface Antigen (HBsAg) Exposure
- Hepatitis C Antibody
- HIV AB Oral Mucosal Transudate
- HIV Ag/Ab Serum
- Influenza
- Influenza Pyrosequencing Assay Clinical
- Lead Filter Paper
- Lead Whole Blood
- Legionella HA
- Lyme Disease EIA
- Measles IgG
- Measles IgM *
- Mumps IgG
- Mumps IgM *
- Mumps PCR *
- Neisseria gonorrhoeae Isolation
- Neisseria Referred Culture
- Norovirus PCR *
- Parasitology
- Pertussis PCR
- Rabies
- Rabies AB Serology
- Rubella IgG

- Rubella IgM *
- Salmonella Serotyping Non-Human *
- Salmonella/Shigella Serotyping Human
- Student/Employee Immune Status Panel
- Syphilis (USR)
- Syphilis DFA
- Syphilis TP-PA * (Must be ordered with Syphilis (USR))
- Syphilis VDRL
- Toxic Shock Testing *
- Trichomonas vaginalis Non-Culture (Fee-for-service only)
- Varicella Zoster IgG
- Viral Culture
- Viral Respiratory PCR Panel

* Prior Approval Required

Though other tests may be listed in the navigation bar, these other tests are being piloted and are not yet ready for use. Additional tests will be brought online soon.

Please note, at this time, ETOR is available only for specimens submitted to the Lansing laboratory.

Where to Get Help

If you need help while using ETOR, please contact the ETOR Application Help Desk:

By email: lims_help@michigan.gov

<u>Notes</u>

<u>Notes</u>