

Orders Status EMR Report

Sample Report – Orders Status Report

The following EMR Report is generated via the EMR reporting tool. It includes all reports that do not have a “completed” status in the Orders table. Tests must be ordered and populate the Orders_ table in the KBM to be included in this report.

order_create_timestamp	Lat Name	Mid Name	Fst Name	DOB	order_actStatus	order_actCode	order_actText
2008-11-06 10:55:53.613	Act		Catherine	19360419	ordered		GLUCOSE
2008-11-06 10:56:17.893	Act		Catherine	19360419	result received		CBC wth Diff
2008-11-06 10:56:18.033	Act		Catherine	19360419	ordered		CBC
2008-11-06 10:56:57.300	Act		Catherine	19360419	ordered	88164	P AP smear
2008-11-25 15:51:40.657	Act		Catherine	19360419	scheduled		
2008-11-25 15:52:05.170	Act		Catherine	19360419			
2008-11-25 16:02:45.077	Act		Catherine	19360419	ordered		CBC
2008-12-15 23:04:24.937	Act		Catherine	19360419	scheduled		Office visit
2009-02-24 14:06:23.573	HENCHAR	J	Pamela	19541213	ordered	74160	CT scan, abdomen, w/ contrast
2009-02-24 14:07:33.263	HENCHAR	J	Pamela	19541213	ordered	88142	P AP, thin prep
2008-10-31 13:43:19.213	Lang		Lana	19830711	ordered		CBC wth Diff
2008-10-31 13:43:19.277	Lang		Lana	19830711	ordered		TOTAL CHOLESTEROL
2008-10-31 13:43:19.430	Lang		Lana	19830711	ordered		TOTAL PROTEIN
2008-10-31 13:43:19.493	Lang		Lana	19830711	ordered		UA w/ MICRO
2009-01-07 13:14:06.190	Test	Girl	Child	19900615	ordered	205	BHCT
2009-01-07 13:19:11.723	Test	Girl	Child	19900615	ordered	85032	CBC, MANUAL
2009-01-07 13:19:11.757	Test	Girl	Child	19900615	ordered	85025	CBC, W/DIFF
2009-01-07 13:19:11.803	Test	Girl	Child	19900615	ordered	87070	CULTURE THROAT
2009-01-21 17:52:48.963	Test		Patient	19710211	ordered	95115	Allergen immunotherapy, one injection

The status of a patient's lab tests are shown in the Lab *Orders* dialog box in NextGen EMR.

The NextGen statuses that display at different times in the first half of the column, before the slash are listed below:

- Pending
- Ordered
- Canceled (displays only when an order is canceled through the Lab Interface)
- Assigned
- Signed-off
- (null) if the order was automatically generated by the Rosetta interface

Report Setup

Go to NextGen EMR -

File > Reports>Generate a Report>By Practice

Templates Tab -

In the Search Template Name field type Order. Select the Order checkbox

In the Fields/Select Field choose the following field boxes:

order_create_timestamp


order_actCode

order_actDescription

order_actStatus


order_actText

order_app.Date



Columns Tab -

Select the items noted below, the order of display can be changed via blue arrows. (Note: Template fields are added to the Columns list for selection only after they are selected in the Template Tab.)




Filter Tab –

Change the Filter Tab as desired. The example below indicates:

The orders created timestamp to occur “before 2/28/2009” other options are:

- On
- Not On
- Before
- After
- On or Before
- On or Between
- Between



(Note: Template fields are added to the Filter Tab for selection only after they are selected in the Template Tab.)