

WIC HemoCue QC Log

Analyzer Type : Hemocue 201+

Site: _____

Control Manufacturer: _____

QC Data (reagent specifics): This area must be completed during the first clinic of each week. If a bottle is opened with a new lot number, a new sheet must also be completed. The package insert will specify the open vial expiration date for high and low controls. Attach a copy of the appropriate package insert to this log.

	Lot number	Close vial expiration date	Open vial expiration date	Expected Range
Low Control				to
High Control				to
Cuvette				na

Q.C. Results

Date	Weekly Cleaning	Self Test (Pass/Fail) (daily)	Low Control (weekly)	High Control (weekly)	Weekly Controls (Pass/Fail)	Initials

Date _____ **Corrective Action:** (use reverse side as needed)

_____/_____/_____
Site Coordinator Review Date

*****To be submitted to Site Coordinator on a monthly basis for review*****