

Resident Name _____

EQUIPMENT AND DEVICE SAFETY LOG – Complete upon initial use and at least every 60 days after the first 5 days following initial use or installation.

The gap between the slats of the bed rail and at the triangular space at the head of the bed can be no greater than 4.5 inches or _____ for this resident.

*The gap between the mattress and the bed rail and between the mattress and the headboard or footboard can be no greater than 2.5 inches or _____ for this resident. These measurements are taken with the resident in bed. Gaps greater those specified require **immediate** corrective action.*

Date	Bed Number	Mattress Number	Device Type	Space at Headboard	Space Beside Right Device	Space Beside Left Device	Space at Footboard	Mattress Affixed	Mechanically Sound Rail	Rail Firmly Attached	Mattress Condition	Bed Frame Condition	Bed Rail Condition	Latch Condition	Corrective Action	Signature