

4.0 Breastfeeding

Effective Date: 2/1/10

4.05A Issuance of Multiple User Breast Pumps

Description of Multi-User Breast Pumps

Name of Pump and Company	Description	Uses
Lactina by Medela	Hospital grade breast pump Allows mother to control suction strength Double or single pumping Multiple-user: each mother needs own kit	Applies to Lactina and Elite: <ul style="list-style-type: none"> • Long-term pumping for NICU • Special needs babies • Mother-baby separation • Breastfeeding problems • Baby not latching • Increasing milk supply • Workplace lactation rooms
Elite by Evenflo/Ameda	Hospital grade breast pump Allows mother to control and adjust suction strength and cycling speed Double or single pumping Multiple-user: each mother needs her own kit	
Nature III by Bailey Medical	Allows mother to control and adjust suction strength Mother manually controls cycling speed Double or single pumping Multiple-user: each mother needs her own kit Overflow filter needs to be changed every 6 weeks	<ul style="list-style-type: none"> • Mother-baby separation • Breastfeeding problems • Baby not latching • Increasing milk supply • Returning to work/school • Pumping multiple times a day

Description of Pedal Pumps

Name of Pump and Company	Manufacturer’s Description	Uses
Pedal Pump by Medela	Medela's Pedal Pump Breastpump uses the strong muscles of the leg and foot to make single or double pumping even easier. This pump requires no electricity and features a unique spring in the pump cylinder, which helps you power the pump. The spring cylinder is attached to an attractive, natural wood pedal and base.	Use in situations where electricity is not available Mother/ infant separation (i.e. return to work/school) Increasing milk supply after breastfeeding is established Mother transitioning from multi-user electric pump Long term use of a double pumping system