

WIC HEMOCUE QUALITY CONTROL (QC) LOG

Michigan Department of Health and Human Services

Site	Analyzer Type
Month	Control Manufacturer

QC Data: Quality Controls must be completed on a weekly basis. An extra line for cuvettes is available if a 2nd bottle is opened during the month. If a new microcuvette bottle with a new lot number is opened, controls must be completed prior to client testing. The package insert for liquid controls will specify open and closed vial expiration dates. Attach a copy of the liquid control insert to this log.

Liquid Controls	Lot Number	Date Opened	Closed Vial Expiration Date	Open Vial Expiration Date	Expected Range
Low Control					
High Control					

Microcuvettes	Lot Number	Date Opened	Expiration Date
Cuvette #1			
Cuvette #2			

QC Results

Date	Low Control (Weekly)	High Control (Weekly)	Weekly Controls (Pass/Fail)	Initials	Cuvette #1 or #2

System Maintenance

(Record error codes below. Use reverse side if additional space is needed.)

Date	Problem / Error Code	Corrective Action	Initials

Site Coordinator Monthly Review _____

Date _____